

Product Name :
Nauga Lab Forced Convection Laboratory Ovens

Product Code :
NLAB-LABORATORYEM40002



Description :

Nauga Lab Forced Convection Laboratory Ovens

Technical Specification :

Forced Convection Laboratory Ovens

Safe, Superior Protection for Sample, User and the Environment.

Multiple redundant over-temperature protection systems to guarantee maximum sample and user protection.

Electronic over-temperature protection built into the microprocessor.

Redundant mechanical over-temperature protection, adjustable from the front, independent from the microprocessor.

Overall temperature protection meets.

Red LED illuminates if external mechanical temperature protection is engaged.

Controller will control temperature at the over temperature setpoint.

All electrical components UL recognized.

Electrical circuit protection in accordance with UL requirements.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003