

**Product Name :**  
Nauga Lab Forced Convection Laboratory Ovens

**Product Code :**  
NLAB-LABORATORYEM40002



**Description :**

Nauga Lab Forced Convection Laboratory Ovens

**Technical Specification :**

Forced Convection Laboratory Ovens

Safe, Superior Protection for Sample, User and the Environment.

Multiple redundant over-temperature protection systems to guarantee maximum sample and user protection.

Electronic over-temperature protection built into the microprocessor.

Redundant mechanical over-temperature protection, adjustable from the front, independent from the microprocessor.

Overall temperature protection meets.

Red LED illuminates if external mechanical temperature protection is engaged.

Controller will control temperature at the over temperature setpoint.

All electrical components UL recognized.

Electrical circuit protection in accordance with UL requirements.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003