

Product Name :
Nauga Lab Hand Grip Dynamometer

Product Code :
NLAB-LABORATORYEM21026



Description :

Nauga Lab Hand Grip Dynamometer

Technical Specification :

These Grip Dynamometer are synonymous with excellent precision.

Engineered following international industry standards, these Grip Dynamometer are widely acclaimed for their precised engineered, accuracy and optimum performance.

Our Grip Dynamometer are available in various technical specifications as per the choices of clients.

Features:

Resistant to weather

Easy to install

Superb quality

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003