

Product Name :
Nauga Lab Heated Bath Circulator

Product Code :
NLAB-LABORATORYEM46017



Description :

Nauga Lab Heated Bath Circulator

Technical Specification :

Heated Bath Circulator

Heated Baths Circulator designed for controlling above ambient temperatures in external instrumentation or for immersion applications.

A very handy tool for general lab requirements and applications.

Compact space saving concept can be placed under or over the work-bench.

Effective water/liquid heated circulator for tighter temperature controls.

Reduces tap water consumption thus reducing high cost of water and waste water.

Waste water does not contain hazardous substances thus eliminating danger to personnel and environment caused by solvents.

Built in circulator pump for efficient circulating processes. Offers highest level of performance, flexibility and control for most demanding applications.

Industrial grade energy efficient heaters reduces power bills drastically

Low noise level.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003