

Product Name :
Naugra Lab High Temperature Oil Bath

Product Code :
NLAB-LABORATORYEM46022



Description :

Naugra Lab High Temperature Oil Bath

Technical Specification :

High Temperature Oil Bath

High Temperature Oil Bath has a chamber made of copper/stainless steel sheet.

Temperature thermostatically controlled up-to 150 C.

Complete with high speed stirrer with speed controller.

Chamber size:-250 x 250 x 250 mm 300 x 300 x 250 mm 355 x 355 x 250 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003