

Product Name :
Nauga Lab Horizontal Laminar Flow Cabinet

Product Code :
NLAB-LABORATORYEM6602



Description :

Nauga Lab Horizontal Laminar Flow Cabinet

Technical Specification :

Horizontal Laminar Flow Cabinet

Specification:-

Horizontal Laminar Air Flow Bench is used for various works that require maintenance of the cleanliness and purity like plant tissue culture, examination and experimentation of precision elect-ronic parts, medical science, pharmacy, agriculture sector and the like.

Designed so as to meet the requirement of US Federal Standard 209 B providing particle free air to meet class 100 conditions.

Cabinet is fabricated of Wooden board laminated with sunmica or Cold rolled steel duly powder coated OR Stainless Steel.

Working table is made of Stainless Steel Sheet.

Side panels and front door are made out of thick transparent plexi glass.

Air drawn through pre-filter is made to pass through highly effective HEPA (High Efficiency Particular Air) filters having efficiency rating as high as 99.99% with cold DOP and 99.97% with hot DOP, thus retaining all air-borne particles of size 0.3 micron and larger.

Heavy Duty Motor operates with minimum noise level.

Working area is illuminated by fluroescent lighting fitted to the unit.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003