

Product Name :
Nauga Lab Hot Air Universal Oven

Product Code :
NLAB-LABORATORYEM38001



Description :

Nauga Lab Hot Air Universal Oven

Technical Specification :

Hot Air Universal Oven

Standard double wall fabrication.

Inner chamber made of richly anodized aluminium or highly polished stainless steel sheet.

Exterior fabricated out of thick mild steel sheet duly finished in stove enamel off-white paint or paste mat finish colour combinations.

Nichrome wire heating elements -provided on three sides attain quick and uniform heating in range of 50oC to 250oC L-shaped thermometer is built-in type.

Control panel is provided with selector switch of high or low rates of power, thermostat control knob and indicators for mains & thermostat.

Supplied with cord and plug.

The equipment is suitable to operate on 220 V AC 50 Hz single phase.

Air Circulating fan.

Digital Temp. Indicator-Cum- Controller.

Digital temp. Indicator

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003