

**Product Name :**  
Naugra Lab Hot Plate Rectangular

**Product Code :**  
NLAB-LABORATORYEM42001



**Description :**

Naugra Lab Hot Plate Rectangular

**Technical Specification :**

Hot Plate Rectangular

Laboratory hot plates are available in round and rectangular shapes.

Body is fabricated out of thick mild steel sheet duly finished in staving enamel paint.

The hot plate made of cast iron, nicely & Precisely machined & smoothed duly finished in heat resistant black paint is firmly mounted on the body.

Heavy duty heating elements are securely layed under the plate to operate on 220 V AC 50 Hz single phase.

Temp. is controlled by a three-heat rotary switch along with an indicator lamp for the mains and cord & plug.

Optional :-

Digital temp. Indicator-cum- controller.

Size In cm	Size in Inches
30x25	(12"x10")
45x30	(18"x12")

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003