

Product Name :
Naugra Lab Hot Plates

Product Code :
NLAB-LABORATORYEM42002



Description :

Naugra Lab Hot Plates

Technical Specification :

Hot Plates

Laboratory hot plates are available in round and rectangular shapes.

Body is fabricated out of thick mild steel sheet duly finished in staving enamel paint.

The hot plate made of cast iron, nicely & Precisely machined & smoothed duly finished in heat resistant black paint is firmly mounted on the body.

Heavy duty heating elements are securely layer under the plate to operate on 220 V AC 50 Hz single phase.

Temp. is controlled by a three-heat rotary switch along with an indicator lamp for the mains and cord & plug.

Specifications as above with 20 cm.(8") dia. Moulded top.

As above but with two plates of 20 cm (8") dia. having independent pilot lamps & switches

Sonic energy regulator can be provided against routine energy regulator.

Digital temp. Indicator- cum- controller

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003