

**Product Name :**  
Naugra Lab Incubator Bacteriological

**Product Code :**  
NLAB-LABORATORYEM37001



**Description :**

Naugra Lab Incubator Bacteriological

**Technical Specification :**

Temp. range from room to 90°C.

Double door, with glass window, Inner chamber of anodized aluminium/stainless steel sheet.

**Specifications:**

Temperature Range  $\pm 0.5^{\circ}\text{C}$

Temperature Uniformity  $\pm 1.0^{\circ}\text{C}$

Temperature Suitable for 18 -24°C, 20-25°C, 30-35°C, 35-37°C, 40-45°C & 55-60°C

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003