

**Product Name :**  
Nauga Lab Incubators Bacteriological

**Product Code :**  
NLAB-LABORATORYEM41001



**Description :**

Nauga Lab Incubators Bacteriological

**Technical Specification :**

Incubators Bacteriological

Standard double wall fabrication.

Inner chamber made out of richly anodized aluminium or highly polished stainless steel sheet and outer made out of thick mild steel sheet finished in white staving enamel paint or mat finished pastel colour combinations.

The gap is filled with special glass wool for temp. insulation.

Double wall door with double window for observation in the chamber.

Temperature Range: 5oC above ambient to 80oC.

Accessories:

Air Circulating fan.

Imported Thermostat "EGO"/ "JUMO" German with the accuracy of  $\pm 1^{\circ}\text{C}$ .

Digital Temp. Indicator

Chamber size in mms Chamber size in Inches

400 x 400 x 400 (16" x 16" x 16")

900 x 900 x 900 (36" x 36" x 36")

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

## **Naugralabequipments**

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com