

Product Name :
Naugra Lab Induction Coil

Product Code :
NLAB-LABORATORYEM21021



Description :

Naugra Lab Induction Coil

Technical Specification :

Induction Coil / Interrupters/ Duboi's Raymond induction coil:

On a teak wood base with graduated scale.

Primary coil with 300 turns and Secondary with 7500 turns with Helmholtz modification of neef's hammer and heavy duty platinum contacts.

Student Induction Coil :

On a teak wood base with graduated scale.

Primary coil with 300 turns and secondary with 6500 turns.

Interrupter working directly off the iron core.

Complete with heavy duty contacts

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003