

Product Name :
Naugra Lab Infra Red Moisture Balance

Product Code :
NLAB-LABORATORYEM51004



Description :

Naugra Lab Infra Red Moisture Balance

Technical Specification :

Infra Red
Moisture Balance

A basic method for determination of moisture content, measure moisture percentage by weight and drying simultaneously.

Infra Red radiations are reflected by a bulb and its temperature is controlled by a regulator.

The balance scale is divided in terms of moisture content on net weight basis.

Specifications:-

Range	
Accuracy	
Dimensions	
Weight Approx.	
Sample Weight	
220V 50Hz	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003