

Product Name :
Naugra Lab Isolated Organ Baths Digital

Product Code :
NLAB-LABORATORYEM21016



Description :

Naugra Lab Isolated Organ Baths Digital

Technical Specification :

Isolated Organ Baths Digital

Versatile and complete organ baths.

Bath chambers made out of thick Transparent acrylic sheets, chemically welded to make them leak-proof.

Temperature is achieved by 250 Watts heating elements and is microprocessor based digital temperature controlled.

A compact control panel has been provided.

Supplied complete with two sets of glass inner vessels, oxygen tubes, Warming coils, fixit levers and

holders.

This bath is termed as 'Two-unit isolated bath'.

The unit is workable on 220V AC 50Hz single phase.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003