

Product Name :
Naugra Lab Laboratory Stirrer

Product Code :
NLAB-LABORATORYEM29001



Description :

Naugra Lab Laboratory Stirrer

Technical Specification :

Heavy & stable cast iron base duly painted with mild steel chrome plated upright rod.

AC/DC motor with stainless steel stirring shaft and blades, variable speed control unit provided.

1/20 HP motor provided for smooth long-time stirring.

Heating of stirrer and stirring space adjustable.

To work on 220 V AC 50 HZ single phase.

Automatic 0-60 minutes timer

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

