


Product Name :
Nauga Lab Magnetic Stirrer with Hot Plate

Product Code :
NLAB-LABORATORYEM42009



Description :

Nauga Lab Magnetic Stirrer with Hot Plate

Technical Specification :

Magnetic Stirrer with Hot Plate

All the specifications are same but in addition it is with hot plate and heat regulator to control the temperature.

The maximum temp. Range up-to 300°C and power rating of heating element is 450 Watt

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003