

**Product Name :**  
Nauga Lab Magnetic Stirrer

**Product Code :**  
NLAB-LABORATORYEM42010



**Description :**

Nauga Lab Magnetic Stirrer

**Technical Specification :**

Magnetic Stirrer with S.S Top

The body is made of rust proof Aluminium sheet.

The Top Plate is made of thick stainless steel of 304 quality.

It is fitted with speed regulator.

The stirring capacity is up-to 2 ltrs.

Complete with 3 Taflon coated magnetic bars of size  $\frac{1}{2}$ ",  $\frac{3}{4}$ " and 1" to work on 220 Volts 50 Hz A.C

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003