

Product Name :
Naugra Lab Magnetic Stirrer

Product Code :
NLAB-LABORATORYEM42010



Description :

Naugra Lab Magnetic Stirrer

Technical Specification :

Magnetic Stirrer with S.S Top

The body is made of rust proof Aluminium sheet.

The Top Plate is made of thick stainless steel of 304 quality.

It is fitted with speed regulator.

The stirring capacity is up-to 2 ltrs.

Complete with 3 Taflon coated magnetic bars of size ½", ¾" and 1" to work on 220 Volts 50 Hz A.C

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003