

**Product Name :**  
Naugra Lab Manometers

**Product Code :**  
NLAB-LABORATORYEM77008



**NAUGRA**®

**Description :**

Naugra Lab Manometers

**Technical Specification :**

All Manometers are fabricated from true-bore Boro-silicate glass, floats are perfectly lapped to eliminate leakage and are fitted with feather weight Stainless Steel capillary

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003