

Product Name :
Naugra Lab Manometers

Product Code :
NLAB-LABORATORYEM77008

NAUGRA®



Description :

Naugra Lab Manometers

Technical Specification :

All Manometers are fabricated from true-bore Boro-silicate glass, floats are perfectly lapped to eliminate leakage and are fitted with feather weight Stainless Steel capillary

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003