


Product Name :
Nauga Lab Membrane Filter Holder All Glass

Product Code :
NLAB-LABORATORYEM65021



Description :

Nauga Lab Membrane Filter Holder All Glass

Technical Specification :

Membrane Filter Holder All Glass

Borosilicate Glass, with Sintered disc for membrane support and aluminium 'duck' clamp.

Ideal for HPLC solvent filtration, Microbiology, etc.

Membrane Filter Holder - All Glass.

For Membrane diameter 47mm.

Funnel Capacity 300 ml

Flask Capacity 1000 ml.

Membrane Filter Holder - All Glass.

For Membrane diameter 25mm.

Funnel Capacity 35 ml.

Flask Capacity 250 ml

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003