

**Product Name :**  
Nauga Lab Membrane Holder

**Product Code :**  
NLAB-LABORATORYEM65026



**Description :**

Nauga Lab Membrane Holder

**Technical Specification :**

Membrane Holder

**Specifications:-**

Diameter 25 / 47 mm.

Buchner funnel of stainless steel.

Vacuum flask of borosilicate glass

Membrane Filter Holder, With Neoprene Lid and Stopper, C-type Handle-clamp.

Membrane Filter Holder, With SS Lid with vent, C-type Handle-clamp.

Membrane Filter Holder, With Neoprene Lid and Stopper and SS Duck Clamp.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003