

Product Name :
Naugra Lab Mini Rotary Flask Shaker

Product Code :
NLAB-LABORATORYEM73001



Description :

Naugra Lab Mini Rotary Flask Shaker

Technical Specification :

Mini Rotary Flask Shaker

Features:-

The complete body is mounted on heavy mild steel angle iron frame making it rigid and sturdy for smooth operations.

Available in Mild Steel painted platform with rubber discs to hold flasks.

Fitted with variable speed motor.

Shaking speed can be adjusted through a speed regulator.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003