

Product Name :
Nauga Lab Muscle Twitch Assembly

Product Code :
NLAB-LABORATORYEM77004



Description :

Nauga Lab Muscle Twitch Assembly

Technical Specification :

Muscle Twitch Assembly

Useful for recording muscle contractions on a smoked paper around cylinder of a drum.

The components of assembly are as under:-

E/8 Recording Drum

Drum Paper Clips

Adjustable stand

Crank Mycrograph

Tuning Fork

Duboi's Raymond Induction coil

Duboi's Raymond Key

Simple Electrode

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003