

Product Name :
Naugra Lab Myographic Lever

Product Code :
NLAB-LABORATORYEM77012



Description :

Naugra Lab Myographic Lever

Technical Specification :

Myographic Lever

Myographic Lever is used to assist with after load screw and fixing nuts to fit on any Myograph

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003