

**Product Name :**  
Naugra Lab Oil Bath (High Temp) with Stirrer

**Product Code :**  
NLAB-LABORATORYEM46008



**Description :**

Naugra Lab Oil Bath (High Temp) with Stirrer

**Technical Specification :**

Oil Bath (High Temp) with Stirrer

Double wall construction

Inner chamber made of highly polished stainless steel or copper sheet and outer made of mild steel duly finished in staving enamel.

Temperature range 50oC to 250oC + 2oC is controlled by a thermostat

Supplied complete with stirrer, diffuser shelves, cover etc.

Panel provided with on-off switch, pilot lamps thermostatic control knob etc.

Chamber Size (in mm)
250x250x250
300x300x250
380x380x250

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003