

| | |
|--|---|
| Product Name : Naugra Lab Open Sided X-Block | Product Code : NLAB-LABORATORYEM77018 |
| Description : Naugra Lab Open Sided X-Block | |
| Technical Specification : Open Sided X-Block Machined from solid Brass rod, richly chromium plated very convenient to hold rods of diameter 3 mm to 13 mm | |

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003