

Product Name :
Naugra Lab Organ Jar Bath Digital

Product Code :
NLAB-LABORATORYEM46023



Description :

Naugra Lab Organ Jar Bath Digital

Technical Specification :

Organ Jar Bath Digital

Made of transparent acrylic sheet, chemically welded to make it perfectly leak-proof.

Capacity of bath 750 ml.

Complete with uprights, organ & oxygen tubes, capillary tubing, frontal lever.

Temperature by immersion heater controlled by a microprocessor based digital temperature controller

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments