

Product Name :
Nauga Lab Oven Vacuum

Product Code :
NLAB-LABORATORYEM40001



Description :

Nauga Lab Oven Vacuum

Technical Specification :

Oven Vacuum

Double wall construction.

Inner round chamber made out of S.S., exterior made of mild steel duly finished in white staving enamel/powder coated paint with mat finished color.

Combination. Temp. range 500C to 1300C.

The body capable of withstanding vacuum of 760 mm hg.

Workable on 200V A.C. 50 Hz single phase.

Supplied without Vacuum pump.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003