

Product Name :
Naugra Lab Pipette Washer

Product Code :
NLAB-LABORATORYEM49001



Description :

Naugra Lab Pipette Washer

Technical Specification :

Pipette Washer

It constitutes of a washing cylinder made of stainless steel sheet highly polished measuring 63 cm.

In diameter and is mounted on a base of 26 cm. diameter.

It is provided with a removable pipette holder.

It can accommodate pipette of 37.5cm length.

It works on the siphoning principle.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003