

**Product Name :**  
Naugra Lab Pipette Washer

**Product Code :**  
NLAB-LABORATORYEM49001



**Description :**

Naugra Lab Pipette Washer

**Technical Specification :**

Pipette Washer

It constitutes of a washing cylinder made of stainless steel sheet highly polished measuring 63 cm.

In diameter and is mounted on a base of 26 cm. diameter.

It is provided with a removable pipette holder.

It can accommodate pipette of 37.5cm length.

It works on the siphoning principle.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003