

Product Name :
Naugra Lab Plasma Thawing Laboratory Bath

Product Code :
NLAB-LABORATORYEM46026



Description :

Naugra Lab Plasma Thawing Laboratory Bath

Technical Specification :

Plasma Thawing Laboratory Bath

Micro Controller based Plasma Thawing Bath, which is a tabletop model, for quick thawing of plasma at 37° C.

In this the temperature is controlled by micro controller based temperature controller.

Our range is designed to safely and reliably thaw fresh frozen plasma (FFP) for the recovery of cryoprecipitate ant hemophilic factor (AHF) and in this the temperature is preset at 4° C and maintained with $\pm 0.2^\circ$ C stability.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003