

Product Name :
Naugra Lab Pneumograph Stethograph

Product Code :
NLAB-LABORATORYEM21010



Description :

Naugra Lab Pneumograph Stethograph

Technical Specification :

Pneumograph Stethograph:

Corrugated rubber tube with side clips, open link chain.

Suitable for use with any tambour

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003