

Product Name :
Nauga Lab Portable Autoclave

Product Code :
NLAB-LABORATORYEM47005



Description :

Nauga Lab Portable Autoclave

Technical Specification :

Portable Autoclave

Portable Autoclave We offer Portable Autoclave which is commonly used in various labs for complete and effective sterilization at the lowest possible cost.

These are known perfect finish, exact in dimensions and sturdy construction.

These are tested well before dispatch and are made from excellent quality material and advanced technology.

These are offered in various shapes and sizes based on the customer requirement.

Features:-

This unit is ever lasting since all the parts are made of non corrosive Stainless Steel Everlasting unit fitted with easily replaceable heater and rubber gasket wire grounded cord and plug Dial Gauge And valve control Complete effective sterilization at the lowest possible coast

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003