

**Product Name :**  
Nauga Lab Propeller Type Mechanical Agitator

**Product Code :**  
NLAB-LABORATORYEM65029



**NAUGRA**®

**Description :**

Nauga Lab Propeller Type Mechanical Agitator

**Technical Specification :**

Propeller Type Mechanical Agitator

Stainless steel rod with propeller type blade.

Fitted with 1/20 HP motor with 4000 RPM with a speed regulator.

Unit fitted on heavy cast-iron base.

It can be raised and lowered in height.

Complete with plug-cord etc.

It works on 220V AC 50Hz single phase.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003