

**Product Name :**  
Naugra Lab Pyrogen Testing Instrument

**Product Code :**  
NLAB-LABORATORYEM75010



**Description :**

Naugra Lab Pyrogen Testing Instrument

**Technical Specification :**

Pyrogen Testing Instrument

Most useful for pyrogen testing of drugs on rabbit, rat's etc. and for measuring the temperature of other animals as great, buffaloes, cow, camels etc.

Microprocessor Based Auto/Manual Scanning LED 4 Digit Display of 8 Channels with Channel Indication, resolution of 0.10C & Accuracy of +1% FS+1 digit.

Measure temperature range from 300C to 450C with the help of Grade "A" PT-100 platinum sensor.

Supplied with 8 nos. of veterinary probes, Operates on 220 V AC+10%, 50 Hz

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003