

Product Name :
Naugra Lab Reciprocating Shaking Machine

Product Code :
NLAB-LABORATORYEM39002



Description :

Naugra Lab Reciprocating Shaking Machine

Technical Specification :

Reciprocating Shaking Machine

Designed for shaking khan racks, bottles and flasks of different sizes by adjusting compartments.

Fitted with heavy duty motor to work on 220 V AC 50 Hz V A single phase.

The rate of motion is Adjustable mechanically by belt tension.

Speed is fixed in same models to approximately 275 or variable as detailed below.

Optional:- 0- 60 minutes timer.

Size of Carrier Stokes/min

20"x13"x4" 250

31"x17"x4" 50-250

36"x18"x4" 50-250

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003