

**Product Name :**  
Naugra Lab Rectangular Muffle Furnace Laboratory Model

**Product Code :**  
NLAB-LABORATORYEM43001



**Description :**

Naugra Lab Rectangular Muffle Furnace Laboratory Model

**Technical Specification :**

Rectangular Muffle Furnace Laboratory Model

Outer shell made out of thick mild steel duly finished in attractive stoving enamel paint.

Fitted with grooved refractory for supporting the heating elements made out of kantar wire and are easily replaceable.

Maximum temp. 1000oC with working temp. 940oC controlled by SUNVIC energy regulator.

Optional :-

Digital temp. Indicator- Cum Controller instead of energy regulator and pyrometer.

**Inner Muffle size**

(in mm)	
---------	--

225x100x100	
-------------	--

250x125x125	
-------------	--

300x125x125	
-------------	--

300'x650'x150	
---------------	--

300'x800'x200	
---------------	--

350'x775'x125	
---------------	--

475'x775'x125	
---------------	--

450'x925'x225	
---------------	--

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003