

Product Name :
Naugra Lab Rectangular Water Bath

Product Code :
NLAB-LABORATORYEM46024



Description :

Naugra Lab Rectangular Water Bath

Technical Specification :

Rectangular Water Bath

Double walled construction. Inner chamber of stainless steel sheet and outer of mild steel sheet duly finished in powder coated.

The gap is filled with glass wool for temp. insulation.

Temperature range 40°C above ambient to 100°C controlled by a thermostat.

Complete with concentric rings, pilot lamp, plug cord etc.

It works on 220V AC 50Hz single phase

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003