

**Product Name :**  
Naugra Lab Rectangular Water Bath

**Product Code :**  
NLAB-LABORATORYEM46024



**Description :**

Naugra Lab Rectangular Water Bath

**Technical Specification :**

Rectangular Water Bath

Double walled construction. Inner chamber of stainless steel sheet and outer of mild steel sheet duly finished in powder coated.

The gap is filled with glass wool for temp. insulation.

Temperature range 40°C above ambient to 100°C controlled by a thermostat.

Complete with concentric rings, pilot lamp, plug cord etc.

It works on 220V AC 50Hz single phase

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003