

Product Name :
Nauga Lab Refrigerated Water Bath

Product Code :
NLAB-LABORATORYEM46020



Description :

Nauga Lab Refrigerated Water Bath

Technical Specification :

Refrigerated Water Bath

Salient Features:-

Versatile usage.

Ergonomic Design

Energy Efficient

Long life

Low Maintenance

Calibration And Protocol Documentation

Features:-

Our system incorporates multiple Devices such as in line process Indicators, environment scanners and temperature or humidity controllers.

Provides Astech cable for direct interface to any dot matrix or laser printer for online or offline data records printing.

It has Bulk data storage capacity with high data retention life.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003