

Product Name :
Nauga Lab Refrigerated Water Bath

Product Code :
NLAB-LABORATORYEM46011



Description :

Nauga Lab Refrigerated Water Bath

Technical Specification :

Refrigerated Water Bath

These Refrigerated Water Baths are fitted with ISI marked CFC free compressor for effective cooling.

Heating is brought about by immersion type heating elements.

Temperature range is from 5 C to 90 C controlled by digital temperature controller with an accuracy of ± 0.5 C .

They are special purpose baths for specific laboratory applications, mounted on a sturdy angle iron frame with inner chamber made of thick S.S and outer of mild steel duly enamel painted.

The gap between the inner and outer wall is filled with high grade glass wool in order to prevent any loss of temperature.

A high-speed stirrer is fitted to maintain temperature uniformity.

The lid of bath is insulated, double walled fixed with hinges and is provided with locking arrangements along with capacity of 28/50/75ltrs (standard)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003