

**Product Name :**  
Nauga Lab Refrigeration Liquid Bath-Cooling Bath

**Product Code :**  
NLAB-LABORATORYEM46009



**Description :**

Nauga Lab Refrigeration Liquid Bath-Cooling Bath

**Technical Specification :**

Refrigeration Liquid Bath-Cooling Bath

Standard double walled construction

Inner chamber made out of 5o C to 90o C and controlled by a sensitive Mercury contact thermometer with an accuracy of + 2oC

Fitted with a stirrer and speed regulator.

Features:-

Corrosion resistance

Stainless finish

Light weighted

(Image approx.)  
Size (in mm)

500x250x250

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003