

**Product Name :**  
Naugra Lab Respiration Pump

**Product Code :**  
NLAB-LABORATORYEM21001



**Description :**

Naugra Lab Respiration Pump

**Technical Specification :**

Single cylinder pump with mechanically operated valve.

All metal construction.

The driving cone has 3 slots of 6", 9", 12" diameter and is accurately balanced.

A 0.25 HP motor workable on 220V AC 50Hz single phase operates with the pump.

The stroke is adjustable and is graduated up to 250 cc.

Motor coupled on adjustable castor wheeled base

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003