

Product Name :
Nauga Lab Rotary Shaker (VDRL Rotator)

Product Code :
NLAB-LABORATORYEM39004



Description :

Nauga Lab Rotary Shaker (VDRL Rotator)

Technical Specification :

Rotary Shaker

Carrier size
13"x 13" fitted with geared motor with concentric motion along with
time switch range 0 to 15 minutes.

Operated on 220/230
Volts.

Available on with timer
& without Timer.

Features:-

Compact design, durable construction

Smooth mixing

Platform with holder clamp and springs

Heavy duty motor

Highly balanced drive mechanism

Non-slip rotating base

Ideal for serological & syphilis tests Uniform speed

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003