

**Product Name :**  
Nauga Lab Rotary Shaker (VDRL Rotator)

**Product Code :**  
NLAB-LABORATORYEM39004



**Description :**

Nauga Lab Rotary Shaker (VDRL Rotator)

**Technical Specification :**

**Rotary Shaker**

Carrier size  
13"x 13" fitted with geared motor with concentric motion along with  
time switch range 0 to 15 minutes.

Operated on 220/230  
Volts.

Available on with timer  
& without Timer.

**Features:-**

Compact design, durable construction

Smooth mixing

Platform with holder clamp and springs

Heavy duty motor

Highly balanced drive mechanism

Non-slip rotating base

Ideal for serological & syphilis tests Uniform speed

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003