

**Product Name :**  
Nauga Lab Rotary Vacuum Evaporator

**Product Code :**  
NLAB-LABORATORYEM71006



**Description :**

Nauga Lab Rotary Vacuum Evaporator

**Technical Specification :**

**Rotary Vacuum Evaporator**

The apparatus consists of an evaporating flask rotated by heavy duty motor in a heating water bath.

Fitted with speed regulator for minor speed adjustments.

Vacuum distillation of solvent at high and low temperature is achieved quickly by use of a heavy duty motor, which prevents foaming and bumping of the liquid (solvent) into the evaporating flask and also distillation under controlled atmosphere.

The evaporating flask is pear shaped to give a large surface area of evaporation.

Evaporating flask is clamped to a glass duct which goes through the drive unit and protrudes about 120mm into the condenser..

Condenser is fitted with standard joint at the top and the bottom for receiver flask.

Supplied complete with Water Bath Digital and fitted with Digital RPM/Speed Meter.

Fitted with Motorized Auto Lift System which operates up or down on touch of a button

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003