

**Product Name :**  
Naugra Lab Skin Fold Caliper

**Product Code :**  
NLAB-LABORATORYEM21027



**Description :**

Naugra Lab Skin Fold Caliper

**Technical Specification :**

The Skinfold Caliper has been used all over the world on projects of international importance and is the most prestigious Caliper available.

It is the only Caliper CE marked under the Medical Devices Directive 93/42/EEC for a Class 1 Device with Measuring Function and is calibrated using masters traceable to National Standards.

Dial Graduation: 0, 20 mm

Measuring Range: 0 mm to 80 mm

Measuring Pressure: 10 gms/mm<sup>2</sup> (constant over range)

Accuracy: 99.00%

Repeatability: 0, 20 mm

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003