

Product Name :
Naugra Lab Slide Warming Table

Product Code :
NLAB-LABORATORYEM42007



Description :

Naugra Lab Slide Warming Table

Technical Specification :

Slide Warming Table

Specification:-

Slide Warming Table are very useful for heating slides to a uniform constant temperature below the melting point of wax.

Top surface made of Stainless Steel (SS-304 grade)

The surface temperature is controlled from ambient to 70 °C.

Slide Warming Table are very useful for heating slides to a uniform constant temperature below the melting point of wax.

Top surface made of Stainless Steel (SS-304 grade)

The surface temperature is controlled from ambient to 70°C.

Supplied complete with pilot lamp, cord and plug but without thermometer.

To work on 220/230 Volts A.C. Supply.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003