

**Product Name :**  
Naugra Lab Slide Warming Table

**Product Code :**  
NLAB-LABORATORYEM42007



**Description :**

Naugra Lab Slide Warming Table

**Technical Specification :**

Slide Warming Table

**Specification:-**

Slide Warming Table are very useful for heating slides to a uniform constant temperature below the melting point of wax.

Top surface made of Stainless Steel (SS-304 grade)

The surface temperature is controlled from ambient to 70 °C.

Slide Warming Table are very useful for heating slides to a uniform constant temperature below the melting point of wax.

Top surface made of Stainless Steel (SS-304 grade)

The surface temperature is controlled from ambient to 70°C.

Supplied complete with pilot lamp, cord and plug but without thermometer.

To work on 220/230 Volts A.C. Supply.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003